


Flame Photometry Rinse Solution Reagecon

[Product in eShop](#)

[Home](#) / [Reagecon_FlamePhotometryStandards](#) / Flame Photometry Rinse Solution Reagecon

Product Data

The Standards are of a very high accuracy, non hazardous and non toxic. Products are manufactured and tested in a GLP (Good Laboratory Practice) environment.

Item No.	Content	Packaging	Price
500.FIRINSE5	500 milliliter	HDPE Bottle	CHF 50.60

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 19.05.2026 at 09:57 o'clock.