

ICP Multi Element Interference Standards Reagecon

[Product in eShop](#)

[Home](#) / [Anion&Cation](#) / ICP Multi Element Interference Standards Reagecon

Product Data

Reagecon have been manufacturing Inorganic Standards, Controls and Calibrators for Spectroscopy for almost two decades. During that time, the company has established itself as the most reliable primary producer of top quality standards. The products are manufactured, tested and stabilised to such a high level, that they can be used on all of the instruments. Reagecon's Multi Element Standards for ICP, ICP-MS are manufactured from ultra-pure metals or salts. Raw materials used for these standards are assayed by titration and ICP-MS prior to manufacture. Separate CRM's are used to control or calibrate the titration and ICP-MS respectively. Products are manufactured gravimetrically using the mass balance approach: 100% less the sum of all impurities (w/w).

Item No.	Content	Packaging	Price
500.REICPINTF10A	100 milliliter	HDPE Triple-leached	CHF 722.73
500.REICPINTF11A	100 milliliter	HDPE Triple-leached	CHF 840.03
500.REICPINTF15A	100 milliliter	HDPE Triple-leached	CHF 941.00
500.REICPINTF15B	100 milliliter	HDPE Triple-leached	CHF 866.54
500.REICPINTF16A	100 milliliter	HDPE Triple-leached	CHF 893.02
500.REICPINTF16B	100 milliliter	HDPE Triple-leached	CHF 1'116.26
500.REICPINTF3A	100 milliliter	HDPE Triple-leached	CHF 336.79
500.REICPINTF4A	100 milliliter	HDPE Triple-leached	CHF 412.46
500.REICPINTF4C	100 milliliter	HDPE Triple-leached	CHF 412.46
500.REICPINTF5A	100 milliliter	HDPE Triple-leached	CHF 723.87
500.REICPINTF8A	100 milliliter	HDPE Triple-leached	CHF 488.15
500.REICPINTF9A	100 milliliter	HDPE Triple-leached	CHF 503.25
500.REICPIS2	100 milliliter	HDPE Triple-leached	CHF 712.60
500.REICPIS5A	100 milliliter	HDPE Triple-leached	CHF 250.85
500.REICPIS7A	100 milliliter	HDPE Triple-leached	CHF 309.15
500.REICPIS7B	100 milliliter	HDPE Triple-leached	CHF 571.90
500.REICPIS7C	100 milliliter	HDPE Triple-leached	CHF 514.64
500.REICPIS7D	100 milliliter	HDPE Triple-leached	CHF 461.63
500.REICPIS7E	100 milliliter	HDPE Triple-leached	CHF 461.63

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 23.04.2026 at 00:41 o'clock.

www.huberlab.ch