

Ion Selective Electrode Standards Reagecon

[Product in eShop](#)

[Home](#) / [Brands](#) / [Reagecon Diagnostics](#) / Ion Selective Electrode Standards Reagecon

Product Data

Ion Selective Electrode (ISE) Standards can be used at suitable concentrations to calibrate or control measurements using Ion Selective Electrodes as the sensing technology. Additionally suitable for Method Validation or Instrument Qualification.

Item No.	Content	Packaging	Price
500.ISEBR5	500 milliliter	HDPE Bottle GL45	CHF 283.35
500.ISECA5	500 milliliter	HDPE Bottle GL45	CHF 118.75
500.ISECL10005	500 milliliter	HDPE Bottle GL45	CHF 233.35
500.ISECL5	500 milliliter	HDPE Bottle GL45	CHF 118.75
500.ISEF0505	500 milliliter	HDPE Bottle GL45	CHF 136.90
500.ISEF08	500 milliliter	HDPE Bottle GL45	CHF 200.70
500.ISEF1	500 milliliter	HDPE Bottle GL45	CHF 291.25
500.ISEF10005	500 milliliter	HDPE Bottle GL45	CHF 258.60
500.ISEF1005	500 milliliter	HDPE Bottle GL45	CHF 223.30
500.ISEF11	1 liter	HDPE Bottle GL45	CHF 364.25
500.ISEF20005	500 milliliter	HDPE Bottle GL45	CHF 200.70
500.ISEF5	500 milliliter	HDPE Bottle GL45	CHF 116.40
500.ISEF50005	500 milliliter	HDPE Bottle GL45	CHF 200.70
500.ISEF5005	500 milliliter	HDPE Bottle GL45	CHF 200.70
500.ISEI5	500 milliliter	HDPE Bottle GL45	CHF 219.60
500.ISEK5	500 milliliter	HDPE Bottle GL45	CHF 219.60
500.ISEMG5	500 milliliter	HDPE Bottle GL45	CHF 194.20
500.ISEN1005	500 milliliter	HDPE Bottle GL45	CHF 187.00
500.ISEN105	500 milliliter	HDPE Bottle GL45	CHF 187.00
500.ISEN5	500 milliliter	HDPE Bottle GL45	CHF 118.75
500.ISENA10005	500 milliliter	HDPE Bottle GL45	CHF 101.65
500.ISENA5	500 milliliter	HDPE Bottle GL45	CHF 219.60
500.ISENH1005	500 milliliter	HDPE Bottle GL45	CHF 323.70
500.ISENH41005	500 milliliter	HDPE Bottle GL45	CHF 113.20
500.ISENH45	500 milliliter	HDPE Bottle GL45	CHF 118.75
500.ISENH5	500 milliliter	HDPE Bottle GL45	CHF 118.75
500.ISEPB5	500 milliliter	HDPE Bottle GL45	CHF 283.35
500.ISES5	500 milliliter	HDPE Bottle GL45	CHF 215.45
500.ISESC5	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 153.20
500.NH3101	100 milliliter	HDPE Bottle	CHF 127.70
500.P10001	1 liter	HDPE Bottle GL45	CHF 302.20

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 23.04.2026 at 00:41 o'clock.

www.huberlab.ch