

Ionic Strength Adjusters Reagecon



Product in eShop



Home / Brands / Reagecon Diagnostics / Ionic Strength Adjusters Reagecon

Product Data

Reagecon's Ionic Strength Adjusters (ISA's) are used in conjunction with Ion Selective Electrodes and Ion Selective Electrode Standards. An ISA, Ionic Strength Adjuster is used to keep the ionic strength of both standard and test high by eliminating matrix effects and suppressing or eliminate interfering ions.

Item No.	Content	Packaging	Price
500.ISACL5	500 milliliter	HDPE Bottle GL45	CHF 118.75
500.ISAI5	500 milliliter	HDPE Bottle GL45	CHF 219.60
500.ISAN5	500 milliliter	HDPE Bottle GL45	CHF 219.60
500.ISANH5	500 milliliter	HDPE Bottle GL45	CHF 377.80
500.ISASC5	500 milliliter	HDPE Bottle GL45	CHF 131.13
500.ISAWH5	500 milliliter	HDPE Bottle GL45	CHF 122.20
500.TISAB1L-500ml	500 milliliter	HDPE Bottle GL45	CHF 33.35
500.TISAB5	5 liter	HDPE Container	CHF 220.60
500.TISAB5L	5 liter	HDPE Jerrycan	CHF 156.30
500.TISABIV	2 liter	Amber Glass Winchester Bottle	CHF 420.00
500.TISAF5	500 milliliter	HDPE Bottle GL45	CHF 116.45
500.TISAF55	5 liter	Bag-in-Box Container	CHF 776.40

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 24.05.2026 at 02:40 o'clock.