

# Indicators Reagecon



Product in eShop

[Home](#) / [Brands](#) / [Reagecon Diagnostics](#) / [Indicators Reagecon](#)

## Product Data

Reagecon's high quality Indicator Solutions are design specifically for manual titrations. High quality indicator solutions such as Reagecon's facilitate end point detection, in a clear and unambiguous manner and are an imperative for accurate manual titrimetry.

**Hazard Statement(s)** H225+H319  
**Class / PG** 3/II  
**UN-No.** UN1993



Item No.	Content	Packaging	Price
500.ADW	2 liter	Amber Glass Winchester Bottle	CHF 156.60
500.ALRED01	100 milliliter	Amber Glass Winchester Bottle GL45	CHF 52.80
500.BRBP00505	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 141.95
500.BRBP0105	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 87.40
500.BRBP0125	125 milliliter	HDPE Bottle	CHF 56.80
500.BRBP0405	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 274.55
500.BRBPIPA0105	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 87.40
500.BRCG0105	500 milliliter	HDPE Bottle GL45	CHF 161.70
500.BRCG05	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 108.95
500.BRCGM05	500 milliliter	Amber Glass Winchester Bottle	CHF 76.75
500.BRP105	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 110.85
500.BRPB040125	100 milliliter	Glass Winchester Bottle	CHF 76.05
500.COR01005	500 milliliter		CHF 76.71
500.COR105	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 89.60
500.CRER0405	500 milliliter		CHF 87.40
500.CVSOLN021	100 milliliter		CHF 95.70
500.DPC05	500 milliliter	Amber Glass Winchester Bottle	CHF 164.15
500.EBB05	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 98.10

500.EBTNA0205	500 gram		CHF 95.70
500.EBTNA0105	500 gram	Amber Glass Bottle	CHF 131.00
500.EOW00051	1 liter	HDPE Bottle GL45	CHF 131.00
500.FEA25	250 milliliter	HDPE Bottle	CHF 162.05
500.FS0101	1 liter	HDPE Bottle GL45	CHF 222.70
500.FS010105	500 milliliter	HDPE Bottle GL45	CHF 113.47
500.FS01015	2 liter	HDPE Container	CHF 146.75
500.FS0102	1 liter	HDPE Bottle GL45	CHF 131.00
500.FS010205	500 milliliter	HDPE Bottle GL45	CHF 109.40
500.IPT01H	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 155.50
500.IPT02H	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 45.95
500.IPT05F	1 liter	Amber Glass Winchester Bottle GL45	CHF 53.55
500.IPT05H	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 49.05
500.IPT05W	2 liter	Amber Glass Winchester Bottle	CHF 112.28
500.IPT10F	1 liter	Amber Glass Winchester Bottle GL45	CHF 105.65
500.IPT10H	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 84.05
500.IPT10W	2 liter	Amber Glass Winchester Bottle	CHF 277.50
500.IPT16W	2 liter	Amber Glass Winchester Bottle	CHF 367.45
500.KR01	100 milliliter	Amber Glass Winchester Bottle	CHF 129.13
500.M004F	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 76.75
500.MGI00505	500 milliliter		CHF 76.75
500.MOS05	500 milliliter	HDPE Bottle GL45	CHF 103.95
500.MTBLU0050250	250 milliliter	Amber Glass Winchester	CHF 54.35
500.MTR050125	125 milliliter	HDPE Bottle	CHF 77.25
500.MTR05025	250 milliliter	HDPE Flip Top Bottle	CHF 103.05
500.MTR060125	125 milliliter	HDPE Bottle	CHF 85.05
500.MTR06025	250 milliliter	HDPE Flip Top Bottle	CHF 278.95
500.OCRSETIN05	500 milliliter	Amber Glass Winchester Bottle	CHF 317.10
500.PCS5	500 milliliter	HDPE Bottle GL45	CHF 73.60
500.PFS1	100 milliliter	Amber Glass Winchester Bottle	CHF 47.85
500.PR505	500 milliliter		CHF 87.40
500.SO0105	500 milliliter	HDPE Bottle GL45	CHF 60.60
500.SO405	500 milliliter		CHF 73.55
500.ST0205	500 milliliter	HDPE Bottle GL45	CHF 74.65
500.ST0505	500 milliliter	HDPE Bottle GL45	CHF 216.90
500.ST1001	1 liter	HDPE Bottle GL45	CHF 98.20
500.ST105	500 milliliter	HDPE Bottle GL45	CHF 314.75
500.TAIND0250	250 milliliter	HDPE Bottle	CHF 67.85
500.TAIND0500	500 milliliter		CHF 74.48
500.TAIND1000	1 liter	HDPE Bottle GL45	CHF 135.80
500.TASHI050H	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 62.70
500.TB04F	500 milliliter	HDPE Bottle GL45	CHF 121.25
500.TBO8F	1 liter	Amber Glass Winchester Bottle	CHF 140.45
500.TECMXI01	100 milliliter	Amber Glass Winchester Bottle	CHF 58.50
500.UN10025	2 liter	Amber Glass Winchester Bottle	CHF 211.05

500.UN1005	50 milliliter	HDPE Bottle	CHF 21.50
500.UN101	100 milliliter	Amber Glass Winchester Bottle	CHF 22.05
500.UN105	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 113.50
BRP0405.0500	500 milliliter		CHF 207.30

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 26.05.2026 at 15:19 o'clock.

www.huberlab.ch