

Combitips® advanced Eppendorf

[Product in eShop](#)

Home / [Labortops_Q1_2026](#) / Combitips® advanced Eppendorf

Product Data

In 1978 Eppendorf revolutionized the industry and laboratory workflow with the introduction of the Combitip/Multipipette system and has been the leader in market innovations for more than 38 years. The Eppendorf Combitips advanced have been completely redesigned and optimized to meet all the needs of any modern laboratory. The Combitips advanced also function according to the positive displacement principle, thus they always dispense the correct volume independently of the density of the liquid and its flow characteristics (e.g. increased vapor pressure or increased viscosity). Working with radioactive or toxic materials also becomes safer due to the hermetically sealed piston preventing aerosol contamination.

- Combitips advanced are the ideal system components for all Eppendorf Multipipettes
- The 9 volume sizes (0.1 mL–?50 mL) offer a maximum range of dispensing volumes with increments as low as 100 nanoliters depending on the Combitip used
- Elongated Combitips in volumes 2.5 mL, 5.0 mL and 10 mL make it possible to reach the bottom of the most common laboratory vessels
- An individual color coding facilitates the quick identification of the desired Combitips
- The unique funnel geometry ensures comfortable handling while preventing damage to gloves
- Combitips advanced are offered in 3 different purity grades to support any application

Packaging: Combitips® advanced and Combitips® advanced PCR Clean are supplied in a dispenser box (100 Combitips® in 4 bags, PCR Clean with zip-lock bag), 25/50 ml Combitips® come in 25 pcs. per box (each containing an adapter).

11.4990.20	0,1 ml	Combitips® advanced Quality	bulk	white	CHF 128.80
11.4990.21	0,2 ml	Combitips® advanced Quality	bulk	light blue	CHF 128.80
11.4990.22	0,5 ml	Combitips® advanced Quality	bulk	purple	CHF 128.80
11.4990.23	1,0 ml	Combitips® advanced Quality	bulk	yellow	CHF 128.80
11.4990.24	2,5 ml	Combitips® advanced Quality	bulk	green	CHF 128.80
11.4990.25	5,0 ml	Combitips® advanced Quality	bulk	blue	CHF 128.80
11.4990.26	10 ml	Combitips® advanced Quality	bulk	orange	CHF 128.80
11.4990.27	25 ml	Combitips® advanced Quality	4 cartons of 25 pieces	red	CHF 151.80
11.4990.28	50 ml	Combitips® advanced Quality	4 cartons of 25 pieces	light grey	CHF 238.56
11.4990.30	- ml	Combitips® advanced assortment pack	contains one Combitips® advanced sample of each size and Adapter advanced 25 and 50 mL	-	CHF 18.85
11.4990.40	0,1 ml	Combitips® advanced Biopur®	individual blister	white	CHF 307.00
11.4990.41	0,2 ml	Combitips® advanced Biopur®	individual blister	light blue	CHF 234.60
11.4990.42	0,5 ml	Combitips® advanced Biopur®	individual blister	purple	CHF 234.60
11.4990.43	1,0 ml	Combitips® advanced Biopur®	individual blister	yellow	CHF 234.60
11.4990.44	2,5 ml	Combitips® advanced Biopur®	individual blister	green	CHF 234.60
11.4990.45	5,0 ml	Combitips® advanced Biopur®	individual blister	blue	CHF 234.60
11.4990.46	10 ml	Combitips® advanced Biopur®	individual blister	orange	CHF 234.60
11.4990.47	25 ml	Combitips® advanced Biopur®	individual blister	red	CHF 277.15
11.4990.48	50 ml	Combitips® advanced Biopur®	individual blister	light grey	CHF 297.85
11.4990.60	0,1 ml	Combitips® advanced PCR Clean	bulk	white	CHF 161.00
11.4990.61	0,2 ml	Combitips® advanced PCR Clean	bulk	light blue	CHF 161.00
11.4990.62	0,5 ml	Combitips® advanced PCR Clean	bulk	purple	CHF 161.00

11.4990.63	1,0 ml	Combitips® advanced PCR Clean	bulk	yellow	CHF 161.00
11.4990.64	2,5 ml	Combitips® advanced PCR Clean	bulk	green	CHF 161.00
11.4990.65	5,0 ml	Combitips® advanced PCR Clean	bulk	blue	CHF 161.00
11.4990.66	10 ml	Combitips® advanced PCR Clean	bulk	orange	CHF 161.00
11.4990.67	25 ml	Combitips® advanced PCR Clean	4 cartons of 25 pieces	red	CHF 190.90
11.4990.68	50 ml	Combitips® advanced PCR Clean	4 carton of 25 pieces	light grey	CHF 205.85

Accessories

Item No.	Description	Price
11.4850.88	Accessories for Combitips advanced® Eppendorf	from CHF 14.15

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 25.05.2026 at 19:40 oclock.

www.huberlab.ch