

## **FILTER-BIO syringe filter, Nylon66, blue, 0.22 $\mu$ m, 13 mm, 100/PK,**



Product in eShop

Home / [FILTER-BIO](#) /

### **Product Data**

#### **Features**

- Color coding: Easier to tell the filter membrane
- Better membrane media: Improved membrane flow rates
- Application Compatibility: Broad range of filtration media meets diverse application needs and best for chromatography application

#### **Application**

- HPLC sample preparation
- Routine QC analysis
- Content uniformity
- Removal of protein precipitate
- Dissolution testing
- Food analysis

- Bioscience analysis
- Environmental samples

Item No.	pore size	diameter	Price
85.3000.00	0.22 µm	13 mm	CHF 73.80
85.3000.01	0.45 µm	13 mm	CHF 73.80
85.3000.02	0.22 µm	25 mm	CHF 92.80
85.3000.03	0.45 µm	25 mm	CHF 92.80
85.3000.04	0.22 µm	30 mm	CHF 142.80
85.3000.05	0.45 µm	30 mm	CHF 142.80

## More product images



\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 01.05.2026 at 01:00 o'clock.