

ToteMax™ Blood Collection Tray



Product in eShop

Home / [Assortment](#) / [Medicine & Healthcare](#) / [Hematology](#) / [Blood Collection](#) / ToteMax™ Blood Collection Tray

Product Data


keep off the floor, work top or bed

Unique phlebotomy tray design and patented hook feature allows the tray to be on or over most bed rail, foot-board and handrails outside a typical patient room.

- Made from food grade resins for a durable non-breakable design
- Deep tray design helps reduce spills
- Dimensions: 432 x 279 x 279 mm

Item No.	colour	Price
15.1021.34	blue	CHF 109.90
15.1023.35	red	CHF 98.08

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 25.04.2026 at 05:06 o'clock.