

Cleanroom gloves KIMTECH™ Sterile G3/G5 Sterile Latex



Product in eShop

Home / [Cleanroom](#) / Cleanroom gloves KIMTECH™ Sterile G3/G5 Sterile Latex

Product Data

Kimtech™ G3 / G5 Sterile Latex Hand Specific Gloves provide a high level of contamination control in critical and controlled cleanroom environments. The gloves are suitable for use in regulated EU GMP cleanrooms of ISO Class 5 or higher, and the sterilised latex ensures safe and seamless protection when and where it counts. A high quality surface with extended texture on the palm and fingertips provides improved grip, tactile sensitivity and comfort in both wet and dry conditions. The natural rubber latex material is sterilised with gamma radiation (at a sterility level (SAL) of 10⁻⁶).

The hand-specific gloves also feature a rolled-edge cuff for added strength and ease of donning, and are double-wrapped to ensure protection and cleanliness. The gloves are also suitable for double donning and single use and are available in a variety of sizes. Rigorous product development and testing processes combined with proven manufacturing processes result in latex cleanroom gloves that fulfil many quality and compliance requirements. The gloves are also supplied with the additional technical information.

- Guaranteed conformity:- PPE Cat. III according to (EU) Regulation 2016/425- EN ISO 374-1:2016 Type C (K) Chemical splash protection- EN 374-4:2003 Resistant to decomposition by chemicals- EN ISO 374-5:2016 Protection against microorganisms and viruses- Sterility assurance level (SAL) 10⁻⁶- Certificate of analysis available online- Certificate of sterility available online
- The fully sterilised latex material offers high contamination protection against proteins, particles, microorganisms and chemical splashes and low endotoxin levels.
- Suitable for use in cleanroom environments according to EU GMP ISO class 5 (level A) or higher. Supplementary technical information is included as proof of compliance.
- The gloves are hand specific, single use and natural coloured with high purity. Trend data is available.
- The texture on the palms and fingertips improves grip and tactile sensitivity for safer and more efficient processes.

The G3 gloves have a higher particle purity than the G5 latex gloves and are also provided with

trend data.

Item No.	size	type no.	Price
3.7519.40	6	G5	CHF 579.99
3.7519.41	6.5	G5	CHF 579.99
3.7519.42	7	G5	CHF 579.99
3.7519.43	7.5	G5	CHF 579.99
3.7519.44	8	G5	CHF 579.99
3.7519.45	8.5	G5	CHF 579.99
3.7519.46	9	G5	CHF 579.99
3.7519.47	10	G5	CHF 569.59
3.7519.50	---	---	CHF 602.50
3.7519.51	---	---	CHF 602.50
3.7519.52	---	---	CHF 602.50
3.7519.53	---	---	CHF 602.50
3.7519.54	---	---	CHF 602.50
3.7519.55	---	---	CHF 602.50
3.7519.56	---	---	CHF 602.50
3.7519.57	---	---	CHF 602.50

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.04.2026 at 23:26 o'clock.