

Gants Kimtech™ G3 Latex

[Product in eShop](#)[Home](#) / [Cleanroom](#) / Gants Kimtech™ G3 Latex

Product Data

Kimtech™ G3 ambidextrous latex gloves offer enhanced tactile sensitivity and performance combined with high contamination control for seamless protection when and where it matters. These natural latex gloves are repeatedly washed in demineralised ultrapure water to ensure consistent control of particulates and extractables. The gloves are then double bagged with polythene materials and a liner in the box to maintain cleanliness. Kimtech™ G3 Latex Cleanroom Gloves are non-sterile and are recommended for ISO Class 3 or higher cleanroom environments. They offer enhanced protection for the wearer and excellent handling in the most demanding applications. A high quality textured surface on the palm and fingertips provides improved grip, tactile sensitivity and comfort in both wet and dry conditions. In addition, the non-powdered design minimises contamination and drying of the wearer's skin. This ensures that cleanroom users are just as well protected as their applications and equipment. The disposable gloves also feature rolled-edge cuffs for easy donning and a better seal with other garments. Rigorous product development and testing processes combined with proven manufacturing processes result in latex cleanroom gloves that fulfil many quality and compliance requirements. In addition, the gloves are supplied with supplementary technical information.

- Guaranteed conformity: PPE category III in accordance with Regulation (EU) 2016/425 and Regulation 2016/425 as amended in the UK, EN ISO 374-1 type C (K) protection against chemical splashes, EN 374-4 resistance to chemical decomposition, EN ISO 374-5, protection against microorganisms and viruses. Certificate of analysis available online
- High-quality latex offers a high level of protection against contamination from microorganisms, chemical splashes, particles and extractable substances
- Suitable for use in ISO Class 3 or higher cleanroom environments Supplementary technical information is included as proof of regulatory compliance
- The gloves are non-sterile, ambidextrous, single-use and natural coloured
- Textured palms and fingertips improve grip and tactile sensitivity for safer and more efficient processes

Item No.	size	Price
3.7521.40	S	CHF 540.75
3.7521.41	M	CHF 540.75

3.7521.42	L	CHF 540.75
3.7521.43	XL	CHF 540.75

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 18.05.2026 at 00:27 o'clock.

www.huberlab.ch