

Stainless steel measuring spoons

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Liquid Handling](#) / [Dose](#) / Stainless steel measuring spoons

Product Data

Stainless steel AISI 304 (1.4301) measuring spoon in a set of 4, as a measuring aid for fast and easy measuring from 1.25 to 15 ml. The volume is engraved on the handle of each measuring spoon.

- Volumes: 1.25 ml, 2.5 ml, 5 ml, 15 ml
- Spoon high-gloss finished, handle matt polished
- Suitable for use with foodstuffs

Item No.	Price
15.5334.80	CHF 19.50

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 27.04.2026 at 13:31 o'clock.