


Replication systems, Nunc™



Product in eShop



Home / [Assortment](#) / [Life Science](#) / [Highthroughput Screening \(HTS\)](#) / [Microplates](#) / Replication systems, Nunc™

Product Data

Designed for the replication of DNA libraries, this system is ideal for transferring small volumes of liquid from plate to plate or to a membrane. Pins transfer 1,0 µl ($\pm 0,1$ µl) of liquid.

- Anodised aluminium replicators with stainless steel pins
- PP copiers
- Dimensions: 121,6 x 83,3 x 135,4 mm

Item No.	execution	Price
7.2503.93	Replicator, 384-pin	CHF 7'087.60
7.2505.20	Replicator, 96-pin	CHF 4'956.55

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 04.05.2026 at 16:30 o'clock.