

## Gas burner - Powerjet WLD-Tec



Product in eShop

[Home](#) / [Assortment](#) / [Research & Development](#) / [Microbiology](#) / [Burners](#) / Gas burner - Powerjet WLD-Tec

### Product Data

Hand-held laboratory gas burner with piezoelectric ignition for liquid gas cartridges. Suitable for flame protection sterilisation in microbiological or biotechnical applications where no central gas supply is available.

- For liquid gas cartridges with 7/16" x 28 UNEF thread, e.g. CG1750
- Ignition of the burner at the touch of a button
- Flame size infinitely adjustable
- Gas supply only when button is pressed
- Locking of the burner button for continuous flaming
- Safe working due to one-handed operation

**Scope of delivery Powerjet:** Laboratory gas burner with gas cartridge adapter (7/16" x 28 UNEF thread). Supplied without gas cartridge. Please order separately

### Specifications Powerjet

- Weight: 300 g
- Gas consumption: 120 g/h
- Flame temperature: 1350 °C
- Warranty: 2 years



Item No.	type no.	Price
13.2053.73	Campingaz gas cartridge CG1750	CHF 94.10
17.3600.11	Powerjet	CHF 460.00

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 12.05.2026 at 09:03 o'clock.

www.huberlab.ch