

Zoom Stereo Microscope SMZ-171 Motic



 [Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Microscopie](#) / [Zoom Stereo Microscope SMZ-171 Motic](#)

Product Data

- Greenough zoom optical system
- Binocular head, 45° inclined, 360° rotating
- Widefield adjustable eyepieces WF10X/23mm
- 6.7:1 Zoom ratio, WD=110mm
- Objective magnification range: 0.75X-5X
- SMZ-171 Pole type stand with large working area and head holder
- Coarse focusing system with tension adjustment
- 3W incident and transmitted LED illumination with independent intensity controls
- Transmitted light with tilting reflector (oblique illumination)
- Option to adapt flexible light guide
- Black/White and frosted glass stage plates, and dust cover

Item No.	description	Price
17.7601.10	Binocular	CHF 2'514.76
17.7601.11	Trinocular	on demand
17.7601.13	---	CHF 1'167.00

[More product images](#)



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 01.05.2026 at 08:19 o'clock.

www.huberlab.ch