

Coliwasa, reusable Wheaton

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Quality Assurance](#) / [Sampling](#) / Coliwasa, reusable Wheaton

Product Data

- Uniquely designed top spin bar and bottom valve plate mechanism provides a positive seal for sample retention along with greater control for releasing the sample
- Manufactured from clear PVC
- Sampling can be achieved up to 20'
- Modular design adds flexibility by enabling the connection of several sections to form a longer sampling unit up to 20' (610cm)
- Provided with a 48" long cleaning brush
- Dimensions (ID x OD): 1 3/8" x 1 5/8" (3.5 x 4.12cm)

Item No.	execution	description	vol.	Price
8.8852.50	Coliwasa Unit	4' (122cm) Coliwasa Unit	760 ml	CHF 725.65
8.8852.52	Coliwasa Unit	6' (183cm) Coliwasa Unit	1100 ml	CHF 971.15
8.8852.53	Extension and Connector	4' (122cm) Extension and Connector	760 ml	CHF 508.45
8.8852.58	Replacement Parts	Bottom Plate and SS Nut	---	CHF 57.55

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.05.2026 at 09:38 o'clock.

www.huberlab.ch