

Lab blenders bag



Product in eShop

[Home](#) / [Assortment](#) / [Medicine & Healthcare](#) / [Sampling](#) / [Sampling bags](#) / [Lab blenders bag](#)

Product Data

- For fast and secure crushing of solid samples
- Special clip for secure fastening of the bags

Item No.	material	l x w	vol.	piece bag/box	information	Price
17.9613.51	PE	180 x 300 mm	400 ml	25/2500	without filter	CHF 720.11
17.9613.52	PE	190 x 300 mm	400 ml	25/2500	full filter	CHF 2'128.23
17.9613.53	PE	190 x 300 mm	400 ml	25/2500	lateral filter	CHF 2'099.34
17.9613.54	PE	380 x 510 mm	3500 ml	100/500	without filter	CHF 318.86

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 04.05.2026 at 11:50 o'clock.

www.huberlab.ch