

## WHEATON - Polypropylene Rack

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Environmental technology](#) / [Sampling & Sample Storage](#) / [WHEATON - Polypropylene Rack](#)

### Product Data

Polypropylene racks for headspace auto-sampler vials.  
Each well has an opening in the bottom to facilitate drainage.

- Autoclavable and stackable
- Chemically resistant
- Include position identifiers

Item No.	number of holes	bore	w x d x h	Price
8.8688.04	48	Ø 15,5 mm	266 x 94 x 28 mm	CHF 307.75
8.8688.05	36	Ø 23,1 mm	322 x 91 x 28 mm	CHF 274.50
8.8688.06	50	Ø 28,1 mm	320x170x35 mm	CHF 445.00
8.8688.10	90	Ø 17,1 mm	320x170x30 mm	CHF 367.50
8.9857.50	96	Ø 9,5 mm	185x125x22 mm	CHF 461.75
8.9858.00	50	Ø 12,5 mm	190 x 100 x 22 mm	CHF 287.75

### More product images



\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 28.04.2026 at 23:49 o'clock.

[www.huberlab.ch](http://www.huberlab.ch)