

## Stand Sartolab® MultiStation Sartorius

[Product in eShop](#)

[Home](#) / [Labortops Q2 2026](#) / [Stand Sartolab® MultiStation Sartorius](#)

### Product Data

Sartolab® MultiStation is a stand specially designed to hold 1 to 6 vacuum filtration units, allowing simultaneous hands-free filtration of up to 6 samples.

- One single vacuum source enables simultaneous filtration of up to 6 samples
- Time-saving (no installation time for each filter unit before use)
- Hands-free filtration

Item No.	type no.	packaging pcs.	Price
12.8900.12	MultiStation	1	CHF 2'590.00 <del>CHF 4'392.00</del>
12.8900.13	Vacuum tube for MultiStation silicone, outer Ø 15 mm, inner Ø 7 mm, length: 200 mm	1	CHF 10.10 CHF 17.15

### More product images



\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

[www.huberlab.ch](http://www.huberlab.ch)