

## Traceable Quartztimer, mit Uhr, mit Standfuss, Clip z. Anstecken und Magnet



Product in eShop

[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Analytics & Measurement](#) / [Time measurement](#) /

### Product Data

- mit LCD-Anzeige
- Speicherfunktion für letzte Programmierung
- Countdown, Alarm und Count-up, Count-down von 24 Std. bis 1 Sek. in 1-Sek.-Schritten
- Genauigkeit 0,01 %
- Gehäuse aus chemikalienresistentem Kunststoff
- Rücken mit Magnet, Clip und Standfuss
- Alarm kann manuell beendet werden oder stellt nach 1 Min. automatisch ab
- mit ersetzbarer 2-Jahresbatterie

| Item No.   | Price     |
|------------|-----------|
| 18.1150.35 | CHF 40.20 |

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 30.04.2026 at 11:00 o'clock.