

Graduated cylinders, SAN, Class B tall shape, with a raised scale Vitlab



Product in eShop

[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Consumables](#) / [Containers](#) / [Cylinders](#) / [Volumetric cylinders](#) / Graduated cylinders, SAN, Class B tall shape, with a raised scale Vitlab

Product Data

Crystal clear.

With molded graduations and ring marks at the primary scale points, calibrated 'In'.

Class B tolerances according to DIN 12681 / ISO 6706.

Hexagonal base with bottom studs provides high stability.

Thermal stress up to 60 °C does not cause tolerance limits to be permanently exceeded.

Suitable for contact with foodstuffs according to regulation (EU) No. 10/2011.

Item No.	division	vol.	h x d	tol. +/-	Price
15.1110.06	1,00 ml	50 ml	199 x 28 mm	1,00 ml	CHF 7.80
15.1110.08	1,00 ml	100 ml	260 x 34 mm	1,00 ml	CHF 11.15
15.1110.10	2,00 ml	250 ml	315 x 47 mm	2,00 ml	CHF 17.65
15.1110.12	5,00 ml	500 ml	350 x 61 mm	5,00 ml	CHF 23.65
15.1110.14	10,00 ml	1000 ml	415 x 76 mm	10,00 ml	CHF 29.50

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.05.2026 at 14:22 oclock.

www.huberlab.ch