

## Multi Timer, 3 Kanal LCD Display 5008CC



 Product in eShop

[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Analytics & Measurement](#) / [Time measurement](#) /

### Product Data

mit LCD-Anzeige, zeigt 3 Count-down oder Alarmzeiten gleichzeitig an, 3 Kanäle können individuell gestartet / programmiert werden

- von 99 Std. 59 Min. 59 Sek. bis 1 Sek. bei 0 ertönt der Alarm, Zeit fertig und die Zeitanzeige beginnt zu zählen
- mit Normalzeit
- Speicherfunktion für die vorangegangene Programmierung
- perfekt für repetitive Tests

Item No.	Price
18.1150.08	CHF 90.00

### Accessories

Item No.	Description	Price
18.1110.39	Batteries for Timer	from CHF 7.10

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 25.04.2026 at 05:06 o'clock.