

HOWORKA - Safety Pipette

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Food technology & cuisine](#) / [Foodtech](#) / [Quality Assurance](#) / [Dilution & Liquid Handling](#) / [HOWORKA - Safety Pipette](#)

Product Data

Out of black rubber/plastic, autoclavable

Laboratory liquid handling device, particularly stable, designed for universal usage. Used for filling and discharging pipettes, burettes and capillary tubes. Also suitable for pipetting and filtering organic solvents and viscous substances.

- The specially designed opening ensures immediate sealing for all commonly used tube diameters
- There are no unnecessary valves, extraneous air cushions or rigid connections

Item No.	suction capacity	execution	Price
15.2176.01	20	small	CHF 40.05
15.2176.02	50	large	CHF 47.41

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 01.05.2026 at 00:58 o'clock.