

Liquid-Sampler single hand operation

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Liquid Handling](#) / [Sampling](#) / [Liquid-Sampler single hand operation](#)

Product Data

Sampler for thin to viscous liquids. Easy use with one-handed operation. With discharge funnel for safe filling of the sample bottles. Material V4A

Type 1 Liquid-Sampler, open with thumb

Type 2 Liquid-Sampler, place onto bottle

Item No.	type no.	total length	volume	for container Ø	Price
15.5326.01	1	460 mm	50 ml	32 mm	CHF 294.60
15.5326.02	1	540 mm	100 ml	32 mm	CHF 324.90
15.5326.03	2	455 mm	50 ml	32 mm	CHF 343.40
15.5326.04	2	535 mm	100 ml	32 mm	CHF 396.65

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 25.05.2026 at 12:55 o'clock.