

Slide boxes

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Medicine & Healthcare](#) / [Microscopy](#) / [Slides](#) / Slide boxes



Product Data

Made of polystyrene, use in freezer down to -80°C and usable up to +70°C. With 50 compartments for standard slides 76 mm x 26 mm x 1 mm

- Archive numbering system in the base
- Stackable
- Colour coded
- Dustproof
- Virtually unbreakable

Item No.	colour	l x w x h	Price
10.6580.00	white	172 x 83 x 30 mm	CHF 6.35
10.6580.01	blue	172 x 83 x 30 mm	CHF 6.35
10.6580.02	green	172 x 83 x 30 mm	CHF 6.35
10.6580.03	red	172 x 83 x 30 mm	CHF 6.35
10.6580.04	yellow	172 x 83 x 30 mm	CHF 6.35
10.6580.05	black	172 x 83 x 30 mm	CHF 6.35
10.6580.10	purple	172 x 83 x 30 mm	CHF 6.35
10.6580.12	orange	172 x 83 x 30 mm	CHF 6.35
10.6580.14	turquoise	172 x 83 x 30 mm	CHF 6.35
10.6599.60	GLW-Divider for 20 slides, up to 1.1 mm width, gray, ZK20	---	CHF 2.35
10.6599.62	GLW-Divider for 25 slides, up to 0.6 mm width, white, ZK25	---	CHF 2.25

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 08.05.2026 at 02:05 o'clock.