

Traceable 4 x stopwatch, with alarm, digital



Product in eShop

[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Analytics & Measurement](#) / [Time measurement](#) /

Product Data

4 Zeiten gleichzeitig

- mit grossem Digitaldisplay
- 4 Kanäle programmierbar von 99 Std. 59 Min. 59 Sek. bis 1 Sek.
- Stoppuhr, Einschaltzeit, Normalzeit
- mit Speicher für vorangegangene Programmierung
- mit magnetischem Rücken - Timer kann an jeder metallischen Fläche befestigt werden
- Genauigkeit 0,01 %
- mit 2 Jahre Silber-Oxyd-Batterie komplett

Item No.	Price
18.1150.04	CHF 88.00

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 25.05.2026 at 12:54 o'clock.