

Ionic Strength Adjusters Reagecon



Produit en ligne

[Home](#) / [Brands](#) / [Reagecon Diagnostics](#) / Ionic Strength Adjusters Reagecon

Données du produit

Reagecon's Ionic Strength Adjusters (ISA's) are used in conjunction with Ion Selective Electrodes and Ion Selective Electrode Standards. An ISA, Ionic Strength Adjuster is used to keep the ionic strength of both standard and test high by eliminating matrix effects and suppressing or eliminate interfering ions.

No. d'article	Contenu	Emballage	Prix
500.ISACL5	500 millilitre	HDPE Bottle GL45	CHF 118.75
500.ISAI5	500 millilitre	HDPE Bottle GL45	CHF 219.60
500.ISAN5	500 millilitre	HDPE Bottle GL45	CHF 219.60
500.ISANH5	500 millilitre	HDPE Bottle GL45	CHF 377.80
500.ISASC5	500 millilitre	HDPE Bottle GL45	CHF 131.13
500.ISAWH5	500 millilitre	HDPE Bottle GL45	CHF 122.20
500.TISAB1L-500ml	500 millilitre	HDPE Bottle GL45	CHF 33.35
500.TISAB5	5 litre	HDPE Container	CHF 220.60
500.TISAB5L	5 litre	HDPE Jerrycan	CHF 156.30
500.TISABIV	2 litre	Amber Glass Winchester Bottle	CHF 420.00
500.TISAF5	500 millilitre	HDPE Bottle GL45	CHF 116.45
500.TISAF55	5 litre	Bag-in-Box Container	CHF 776.40

* Les prix sont sans engagement et s'entendent comme prix d'achat en francs suisses sans la TVA et à l'exclusion de tous autres frais, charges et taxes. Les prix affichés dans l'eShop peuvent être différents du fichier PDF en raison des mises à jour régulières.

** Veuillez noter que lors de la commande de produits chimiques et de détergents, des frais de transport et d'emballage pour les marchandises dangereuses ainsi que des frais légaux sont facturés. Ces frais vous seront indiqués en détail sur la confirmation de commande que vous recevrez en complément de l'accusé de réception de la commande.

*** Vous trouverez de plus amples informations, telles que des informations techniques et des fiches de données de sécurité, dans notre boutique en ligne.

**** Le fichier PDF a été créé sur www.huberlab.ch le 10.05.2026 à 13:52 heures.