

# Spectrophotometry Standards Reagecon



Produit en ligne

[Home](#) / [Brands](#) / [Reagecon Diagnostics](#) / Spectrophotometry Standards Reagecon

## Données du produit

Reagecon's Spectrophotometry Standards are suitable for use on all UV-VIS Spectrophotometers. The range of standards includes: Linearity, Wavelength, Stray Light, Bandwidth and Resolution Standards.

No. d'article	Contenu	Emballage	Prix
500.CAFM1000	1 millilitre	Amber Glass Ampoule	CHF 213.05
500.CAFSETWAT6	30 millilitre	Amber Glass Ampoule	CHF 1'056.00
500.CAFSTD100005	5 millilitre	Amber Glass Ampoule	CHF 184.70
500.CAFSTD10002	2 millilitre	Amber Glass Ampoule	CHF 84.30
500.CAFSTD10005	5 millilitre	Amber Glass Ampoule	CHF 188.35
500.CAFSTD140002	2 millilitre	Amber Glass Ampoule	CHF 103.25
500.CAFSTD1M05	5 millilitre	Amber Glass Ampoule	CHF 443.65
500.CAFSTD200005	5 millilitre	Amber Glass Ampoule	CHF 437.55
500.CAFSTD20005	5 millilitre	Amber Glass Ampoule	CHF 188.35
500.CAFSTD220002	2 millilitre	Amber Glass Ampoule	CHF 103.25
500.CAFSTD25005	5 millilitre	Amber Glass Ampoule	CHF 188.35
500.CAFSTD2K002	2 millilitre	Amber Glass Ampoule	CHF 278.55
500.CAFSTD2M01	100 millilitre	Amber Glass Bottle	CHF 218.05
500.CAFSTD300002	2 millilitre	Amber Glass Ampoule	CHF 103.25
500.CAFSTD40005	5 millilitre	Amber Glass Ampoule	CHF 184.70
500.CAFSTD4M001	1 millilitre	Amber Glass Ampoule	CHF 237.60
500.CAFSTD500005	5 millilitre	Amber Glass Ampoule	CHF 437.55
500.CAFSTD50005	5 millilitre	Amber Glass Ampoule	CHF 341.65
500.CAFSTD5005	5 millilitre	Amber Glass Ampoule	CHF 188.35
500.CAFSTD60002	2 millilitre	Amber Glass Ampoule	CHF 103.25
500.RSPEC-CP015	1 pièce	Quartz Glass Cuvette Permanently Sealed	CHF 795.85
500.RSPEC-CP0151	100 millilitre	Amber Glass Winchester Bottle	CHF 297.55
500.RSPEC-CP036	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 864.95
500.RSPEC-CP0361	100 millilitre	Amber Glass Winchester Bottle	CHF 348.85
500.RSPEC-CP0362	100 millilitre	Amber Glass Winchester Bottle	CHF 188.55
500.RSPEC-CP037	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 864.95
500.RSPEC-CP0371	100 millilitre	Amber Glass Winchester Bottle	CHF 188.55

500.RSPEC-CP0372	100 millilitre	Amber Glass Winchester Bottle	CHF 157.15
500.RSPEC-EP0002	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 1'117.15
500.RSPEC-EP001500	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 1'400.55
500.RSPEC-EP0030	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 951.55
500.RSPEC-EP003600	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 1'566.05
500.RSPEC-EP0058	3 set	Quartz Glass Cuvette Permanently Sealed	CHF 2'854.90
500.RSPEC-EP0059	2 set	Quartz Glass Cuvette Permanently Sealed	CHF 1'834.47
500.RSPEC-EP0062	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 963.65
500.RSPEC-EP00621	100 millilitre	Amber Glass Winchester Bottle	CHF 373.85
500.RSPEC-EP0062X2BU	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 1'027.30
500.RSPEC-EP00631	100 millilitre	Amber Glass Winchester Bottle	CHF 202.10
500.RSPEC-EP0064	1 pièce	Quartz Glass Cuvette Permanently Sealed	CHF 908.20
500.RSPEC-EP006400	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 1'086.10
500.RSPEC-EP00641	100 millilitre	Amber Glass Winchester Bottle	CHF 170.20
500.RSPEC-EP0071	1 pièce	Quartz Glass Cuvette Permanently Sealed	CHF 704.17
500.RSPEC-EP0072	1 pièce	Quartz Glass Cuvette Permanently Sealed	CHF 662.00
500.RSPEC-EP00731	100 millilitre	Amber Glass Winchester Bottle	CHF 202.10
500.RSPEC-EP00741	100 millilitre	Amber Glass Winchester Bottle	CHF 202.10
500.RSPEC-USP01	1 pièce	Quartz Glass Cuvette Permanently Sealed	CHF 1'027.30
500.RSPEC-USP0221	100 millilitre	Amber Glass Winchester Bottle	CHF 116.40
500.RSPEC-USP0261	100 millilitre	Amber Glass Winchester Bottle	CHF 116.40
500.RSPEC-USP0511	100 millilitre	Amber Glass Winchester Bottle	CHF 300.80
500.RSPEC0001	1 pièce	Quartz Glass Cuvette Permanently Sealed	CHF 1'013.90
500.RSPEC00011	100 millilitre	Amber Glass Winchester Bottle	CHF 278.35
500.RSPEC0008	1 pièce	Quartz Glass Cuvette Permanently Sealed	CHF 1'013.90
500.RSPEC00081	100 millilitre	Amber Glass Winchester Bottle	CHF 150.45
500.RSPEC0015	1 pièce	Quartz Glass Cuvette Permanently Sealed	CHF 795.85
500.RSPEC0015060	60 millilitre	HDPE Poly Bottle	CHF 157.15
500.RSPEC00151	100 millilitre	Amber Glass Winchester Bottle	CHF 296.90
500.RSPEC00152	100 millilitre	Amber Glass Winchester Bottle	CHF 192.30
500.RSPEC0020	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 938.21
500.RSPEC00201	100 millilitre	Amber Glass Winchester Bottle	CHF 167.85
500.RSPEC0022	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 938.21

500.RSPEC00221	100 millilitre	Amber Glass Winchester Bottle	CHF 310.54
500.RSPEC0023	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 938.21
500.RSPEC00231	100 millilitre	Amber Glass Winchester Bottle	CHF 167.85
500.RSPEC0024	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 938.21
500.RSPEC0024060	60 millilitre	HDPE Poly Bottle	CHF 137.13
500.RSPEC00241	100 millilitre	Amber Glass Winchester Bottle	CHF 310.54
500.RSPEC0025	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 938.21
500.RSPEC00251	100 millilitre	Amber Glass Winchester Bottle	CHF 167.85
500.RSPEC0026	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 938.21
500.RSPEC00261	100 millilitre	Amber Glass Winchester Bottle	CHF 167.85
500.RSPEC00311	100 millilitre	Amber Glass Winchester Bottle	CHF 192.30
500.RSPEC0036	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 991.10
500.RSPEC00361	100 millilitre	Amber Glass Winchester Bottle	CHF 192.30
500.RSPEC0037	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 991.10
500.RSPEC00371	100 millilitre	Amber Glass Winchester Bottle	CHF 192.30
500.RSPEC0038	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 991.10
500.RSPEC00381	100 millilitre	Amber Glass Winchester Bottle	CHF 355.80
500.RSPEC0039	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 991.10
500.RSPEC00391	100 millilitre	Amber Glass Winchester Bottle	CHF 192.30
500.RSPEC0040	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 991.10
500.RSPEC00401	100 millilitre	Amber Glass Winchester Bottle	CHF 192.30
500.RSPEC0041	1 set	Quartz Glass Cuvette Permanently Sealed	CHF 991.10
500.RSPEC0041060	60 millilitre	HDPE Poly Bottle	CHF 157.15
500.RSPEC00411	100 millilitre	Amber Glass Winchester Bottle	CHF 192.30
500.RSPEC0051	1 pièce	Quartz Glass Cuvette Permanently Sealed	CHF 357.75
500.RSPEC00511	100 millilitre	Amber Glass Winchester Bottle	CHF 355.80
500.RSPEC0052	1 pièce	Quartz Glass Cuvette Permanently Sealed	CHF 357.75
500.RSPEC00521	100 millilitre	Amber Glass Winchester Bottle	CHF 167.85
500.RSPEC00531	100 millilitre	Amber Glass Winchester Bottle	CHF 192.30
500.RSPEC00541	100 millilitre	Amber Glass Winchester Bottle	CHF 192.30
500.UV2801001	100 millilitre	Amber Glass Winchester Bottle	CHF 60.90

\* Les prix sont sans engagement et s'entendent comme prix d'achat en francs suisses sans la TVA et à l'exclusion de tous autres frais, charges et taxes. Les prix affichés dans l'eShop peuvent être différents du fichier PDF en raison des mises à jour régulières.

\*\* Veuillez noter que lors de la commande de produits chimiques et de détergents, des frais de transport et d'emballage pour les marchandises dangereuses ainsi que des frais légaux sont facturés. Ces frais vous seront indiqués en détail sur la confirmation de commande que vous recevrez en complément de l'accusé de réception de la commande.

\*\*\* Vous trouverez de plus amples informations, telles que des informations techniques et des fiches de données de sécurité, dans notre boutique en ligne.

\*\*\*\* Le fichier PDF a été créé sur [www.huberlab.ch](http://www.huberlab.ch) le 08.05.2026 à 22:11 heures.

[www.huberlab.ch](http://www.huberlab.ch)