

Buffer Standard pH 4.00/7.00 at 20/25°C Colour Coded ISO 17034 Certified Reference Material Reagecon



Produit en ligne

Home / [Labortops Q4 2025](#) / Buffer Standard pH 4.00/7.00 at 20/25°C Colour Coded ISO 17034 Certified Reference Material Reagecon

Données du produit

Reagecon's Buffer Standard Certified Reference Material is colour coded for ease of use, produced under our ISO 17034 (A2LA Ref: 6739.01) Accreditation and is NIST traceable.

No. d'article	Contenu	Emballage	Prix
500.CRM1040525C	1 pièce		CHF 23.00 CHF 32.60
500.CRM10405C	1 pièce		CHF 22.50 CHF 31.60
500.CRM1070525C	1 pièce		CHF 23.00 CHF 32.60
500.CRM10705C	1 pièce		CHF 22.50 CHF 31.60

* Les prix sont sans engagement et s'entendent comme prix d'achat en francs suisses sans la TVA et à l'exclusion de tous autres frais, charges et taxes. Les prix affichés dans l'eShop peuvent être différents du fichier PDF en raison des mises à jour régulières.

** Veuillez noter que lors de la commande de produits chimiques et de détergents, des frais de transport et d'emballage pour les marchandises dangereuses ainsi que des frais légaux sont facturés. Ces frais vous seront indiqués en détail sur la confirmation de commande que vous recevrez en complément de l'accusé de réception de la commande.

*** Vous trouverez de plus amples informations, telles que des informations techniques et des fiches de données de sécurité, dans notre boutique en ligne.

**** Le fichier PDF a été créé sur www.huberlab.ch le 19.05.2026 à 09:54 heures.