

Machine à emballer sous vide V.200



Produit en ligne



Home / [Assortiment](#) / [Alimentation & Cuisine](#) / Machine à emballer sous vide V.200

Données du produit

Functions

- Fully automatic: Automatic capping vacuuming and welding process. Monitored by sensor control
- Professional vacuuming and welding of food, industrial parts and other items
- Bag or vacuum sealing and welding rollers, reversible at any time on the container mode
- Individually adjustable sealing time - so you can use foils up to 200 (industrial strength) with a width of up to 350 mm
- Manufacture of foil bags of individual length (continuous rolls can be found in 4 widths in the accessories section)
- Soft vacuum: possible Precise pressure regulation (additional) for vacuum sealing of pressure-sensitive foods
- No pressure on regulatory stop-pulse function or manual welding as part of other devices in this sector
- Vacuuming of all vacuum containers, pots, jars and twist-off jars and containers by suction
- Welding of plastic bags up to 200 to be used as a pure ideal sealer
- Vacuum packing and sealing of liquids (with the L+ pressure regulation or pre-frozen)
- Vacuum packing of smooth, standard vacuum-sealed foil available from 120 thickness (cost savings)!

Vacuum
processes
in succession

No. d'article	Prix
18.0301.03	CHF 557.43

Accessoires

No. d'article	Description	Prix
18.0301.12	Rouleaux continus d'emballage sous vide E-Vac	dès CHF 36.64
18.0301.50	Sacs sous vides G-Vac	dès CHF 30.15
18.0301.53	Sacs sous vides RS-Vac	dès CHF 25.70
18.0301.20	Sacs sous vides R-Vac	dès CHF 21.32
18.0301.70	Récipients Sous vide ronde New-line	dès CHF 20.32
18.0301.57	Sacs sous vides i-Vac	dès CHF 155.04
18.0301.90	Récipients sous vide en acier inoxydable ES-line	dès CHF 57.86
18.0301.80	Récipients Sous vide carré New-line	dès CHF 28.26
18.0301.93	Couvercles universal sous vide Top	dès CHF 12.20
18.0302.04	Pompe manuelle Easy-Pump	dès CHF 17.47
18.0302.03	Utensile d'aspiration pour couvercles et récipients	dès CHF 24.12
18.0302.09	Régulateur de dépressurisation L+	dès CHF 15.81
18.0302.02	Bouchons de bouteille sous vide	dès CHF 14.73
18.0302.12	Récipient sous vide aux normes gastronomiques	dès CHF 178.98

Plus de photos de produits



* Les prix sont sans engagement et s'entendent comme prix d'achat en francs suisses sans la TVA et à l'exclusion de tous autres frais, charges et taxes. Les prix affichés dans l'eShop peuvent être différents du fichier PDF en raison des mises à jour régulières.

** Veuillez noter que lors de la commande de produits chimiques et de détergents, des frais de transport et d'emballage pour les marchandises dangereuses ainsi que des frais légaux sont facturés. Ces frais vous seront indiqués en détail sur la confirmation de commande que vous recevrez en complément de l'accusé de réception de la commande.

*** Vous trouverez de plus amples informations, telles que des informations techniques et des fiches de données de sécurité, dans notre boutique en ligne.

**** Le fichier PDF a été créé sur www.huberlab.ch le 10.05.2026 à 12:37 heures.